



*numo
arms*

User Manual

**Numo
Arms**

**INNOVATIVE
SOLUTIONS**

KCare
Healthcare Solutions

Safety Precautions

The **KA410NAS Numo Toilet Arms** are designed and manufactured with safety and function in mind for the user.

Ensure that the Numo toilet arms are securely attached to the toilet to ensure a stable platform for the user.

Users **MUST** read and follow the User Manual to ensure that all instructions are followed and not compromised.


Specifications

KA410NAS

Overall Width	625mm
Overall Depth	385mm
Overall Height	250mm – 280mm
Seat Width	480mm
Product weight:	6kg
Max User Weight	150kg

Materials


Mild steel base - white nylon plastic coated. Nylon reinforced plastic armrest




Assembly

The KA410NAS Numo toilet arms require assembly prior to use.

1. Slide the plastic arms over the base plate upstands. The arms may have to be tapped down into place to a position where the bottom screw hole aligns with the bottom hole of the base plate upstand.




2. Insert the bolts making sure the countersunk heads are located on the inside of the arms and tighten nut over washer. **Do not over tighten!**



The **KA410NAS Numo toilet arms** is now ready for installation

Installation

Remove the existing toilet seat and place the Numo Arms base plate over the pan, using the two rear holes to align with the pan holes. Locate the seat back over the base plate and pan securing the seat in the normal way with secure fixing bolts holding the seat.



NOTE:

Some older style seats may have a skirt to the sides which may prevent the seat sitting down flat on the front of the pan with the Numo Toilet Arm installation. You may choose to purchase a new seat (making sure it is flat on the underside) from your plumbing supplier or hardware store.

Usage

The **KA410NAS Numo Arms** should be installed and used as per manufacturers recommendations to ensure safe and long lasting trouble free usage of the product.

Before Use

Inspect and ensure the Numo Arms are securely attached to the toilet with the base plate being in even contact with the pan.

Attention!

Heavy or uneven weight distribution may cause stress to the nylon or metal fixing bolts and may eventually result in damage to the product.

Using the Numo Arms Correctly

The weight capacity of this product may be reduced if the user applies their weight unevenly, forcefully or positions the product unevenly.



Correct



Wrong




Warning: Uneven or unsecure fixing of the Numo Arms to the toilet base may result in injury to the user or damage to the Numo arms or toilet.

Inspection

The **KA410NAS Numo Arms** should be inspected on a daily basis for damage.

Check The following:

- Securely fixed to toilet
- Arms
- Base Plate
- Fasteners



Warning: If any part of the Numo Arms is damaged or unsecure, it should be removed from service and not be used until repaired or replaced.

Cleaning

The **KA410NAS Numo Arms** should be washed weekly with a good quality non-abrasive bathroom cleaner. Periodically remove the seat and numo toilet arms and clean between the base plate and pan.

- Do not use harsh chemicals or abrasive cleaning products.
- Do not high pressure clean the numo toilet arms.
- Do not autoclave.

Servicing

We recommend a 6-month preventative maintenance service on the **KA410NAS Numo Arms**.

Please contact K Care Healthcare Solutions service department or your local dealer to arrange service.

Preventative maintenance service includes:


- Examine frame for structural problems
- Examine plastic arms for damage and wear
- Examine nuts and bolts for damage and wear

Numo Arms are to be inspected and repaired if necessary. (see next page for spare parts information)

Important Note: Only trained/qualified technicians should service numo toilet arms.

Spare parts identification

Please refer to the below diagram to identify spare parts



Warranty Against Defects

12 Month Product Warranty Period Applies

1. K Care Healthcare Solutions provides consumers with the following warranty in relation to this Product, in addition to complying with the requirements of any relevant legislation, including the Competition and Consumer Act 2010 (Cth) in Australia and the Consumer Guarantees Act 1993 in New Zealand (the Acts), except where a New Zealand consumer acquires the Product for the purposes of a business.
 - (a) Form means a warranty claim form issued by K Care Healthcare Solutions in respect of Products.
 - (b) K Care Healthcare Solutions or our means K Care Healthcare Solutions, ABN 47 159 431 099 836A Fifteenth St Mildura Vic 3500.
 - (c) Products means the goods manufactured by K Care Healthcare Solutions (including products manufactured by its contract manufacturers) set out in the table above;
 - (d) Material means a material or component used by K Care Healthcare Solutions in the manufacture of the Products;
 - (e) Retailer means the authorised dealer of Products from whom the Product was purchased; (f)Warranty Period means the warranty period set out in the table above,commencing from the date of purchase of the relevant Product. If the Product or part (as the case may be) is repaired or replaced, there will be no extension to the original warranty period; and
 - (g) Workmanship means the handling, assembly and manufacturing processes performed by or on behalf of K Care Healthcare Solutions in order to manufacture the Products.
2. In this warranty, we have used the following definitions:
 3. K Care Healthcare Solutions warrants that for the duration of the Warranty Period, all Products will be free of faults arising from defects in Workmanship or Materials, on the terms and conditions set out in this warranty.
 4. K Care Healthcare Solutions undertakes that if during the Warranty Period any Product, or any part of a Product, has failed to operate correctly due to faulty Workmanship or defective Material, it will repair or replace the Product or part (as the case may be) free of charge (excluding freight fees if applicable) provided that the following procedure is met:

Improving your life everyday™

At Your Service

Sales, administration
and general enquiries

National 1300 783 783

Your Customer Service Team
customerservice@kcare.com.au

kcare.com.au

All information is considered to be true and correct at the date of publication. Changes in circumstances after the time of publication may impact on the accuracy of the information.

